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
thermoscientific

TECHNICAL NOTE 73768

Quantification of free metanephrines in human plasma by LC-MS/MS for clinical research

Authors: Mariana Baroenas and Marie Calvet
Thermo Fisher Scientific, Les Ulis, France

Keywords: Metanephrines, online analysis, plasma, mass spectrometry, TSO Altis

Application benefits

- Accurate and confident results obtained by implementation of a comprehensive kit for sample preparation
- Robust, sensitive LC and MS platforms enable increased confidence in data

Goal
Implementation of an analytical method for the quantification of free metanephrines in human plasma on a Thermo Scientific™ TSO Altis™ triple quadrupole mass spectrometer.

Introduction
Plasma free metanephrines (PFM) are the most specific indicators for the diagnosis of extra-adrenal chromaffin tumors. Because catecholamines are metabolized within chromaffin cells to metanephrine or normetanephrine, these metabolites can be used to diagnose pheochromocytoma. The measurement of PFM (metanephrine (M), normetanephrine (NM), and 3-methoxytyramine (MT)) is challenging because of their polar nature, their low molecular weight, and their low physiological concentration in human plasma.

In this report, an analytical method for clinical research for the quantification of metanephrines in human plasma is reported. Samples were pre-processed by protein precipitation followed by online solid-phase extraction (SPE) and chromatographic separation on a Thermo Scientific™ Vanquish™ Flex Binary UHPLC system. Detection was performed on a TSO Altis triple quadrupole mass spectrometer with heated electrospray ionization (HESI) operated in positive ion mode by selected reaction monitoring (SRM). Method performance was evaluated using the ClinMass™ LC-MS/MS Complete Kit for Free Metanephrines in Plasma – online analysis from RECIPE Chemicals + Instruments GmbH (Munich, Germany) in terms of linearity of response within the calibration ranges, accuracy, and intra- and inter-assay precision for all analytes.

ThermoFisher SCIENTIFIC

PAL SYSTEM
Ingenious sample handling

LC/MS Application Note

LIFE SCIENCE **CLINICAL**


Automated Sample Preparation for Therapeutic Drug Monitoring by LC-MS

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Park Station | Peltier Stack 2DW (4°C) | Fast Wash Station | Tray Holder | Solvent Module | Peltier Stack 2DW (37°C) | Valve Drive

Work flow

Sample pretreatment:

150 μ l Mixture IS/P*

Precipitation:

(A) 50 μ l Sample (calibrator, control, patient)
Add with shaking (1000 rpm)!

shake (10 min, 1000 rpm)

Transfer:

(A) place 96-well filter plate A on 96-well plate B

(B)

centrifuge (10 min, 2000 x g)

LC-MS/MS Analysis:

(B) remove plate A, cover plate B (protective sheet)

Inject 10–50 μ l** of the filtrate

+RECIPE

+RECIPE

ClinMass®
TDM 200 Kit System

Therapeutic Drug Monitoring (TDM)
by LC-MS/MS

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